


# SYRINGE PUMP CONTROLLER

## SXPC Series



### Features

- This controller can be used for both, SBP and SCP Series.
- By connecting the controller to a PC via USB, the operating speed and volume of the pump can be stored in memory and the pump can be controlled manually or from a PC.
- To use the SXP-APS software for your computer, please download it from the QR code on the right.



### Specifications

Item	Specifications
Input Voltage	DC24 V
Power Consumption	6 W
Compatible pumps ※	SBP-100G Series / SCP-1000G Series
Number of pumps that can be connected	1 Unit
Dimensions	W125×D85×H26.8 mm
Weight	150 g
Compatible OS	Windows® 10
Accessories	AC adapter

※Compatible Models	Operating speed	Travel Length
SBP-100G Series	0.01 ~ 20.00 μL/s	0.01 ~ 120.00 μL
SCP-1000G Series	0.10 ~ 500.00 μL/s	0.10 ~ 1200.00 μL

※The official name of Windows® is Microsoft® Windows® Operating System. Windows is a registered trademark or trademark of Microsoft Corporation in the United States and other countries.

Note: Details including specifications may change at any time without notification.

### TFS (Takasago Fluidic Systems)

1900 West Park Drive, Suite 280, Westborough, MA 01581 U.S.A.

Tel: +1-508-983-1434 FAX: +1-508-983-1401

E-mail: info@takasago-fluidics.com URL: https://www.takasago-fluidics.com/